LTV[®] 950 Ventilator

Light Enough to Move Well Advanced Enough to Breathe Well

The LTV® 950 ventilator weighs only 13.4 pounds yet provides many advanced functions required for today's therapy needs – pressure control, pressure support and flow-triggering. This sophisticated technology combined with enhanced portability enables the user to be mobile and comfortable.







LTV[®] 950 Ventilator

Specifications

Physical Specifications:

Weight: Height Width: Depth:

13.4 lbs. (6.1 kg) 3 inches (8 cm) 10 inches (23 cm) 12 inches (30 cm)

Controls:

Power: Ventilation Modes: Breath Types:

Breath Rate: Tidal Volume: Inspiration Time: PC/PS/Spontaneous Flow: **Pressure Control: Pressure Support:** Sensitivity: **Bias Flow** PEEP/CPAP: Control Lock: Manual Breath:

Variable Alarms:

Apnea Interval: High Pressure Limit: Low Peak Pressure: Low Minute Volume: Alarm Volume: Alarm Silence/Reset: Ventilator Inoperative: High Rate Alarm: High PEEP Alarm:

Fixed Alarms:

Disconnect/Sense Line Low & Lost External Power Low & Empty Internal Battery High & Low Oxygen Inlet Pressure

Monitors and Indicators:

Peak Inspiratory Pressure: Mean Airway Pressure: PEEP: Total Breath Rate: Airway Pressure Display: **Exhaled Tidal Volume:** Total Minute Volume: IE Ratio: Calculated Peak Flow: Patient Effort:

14.45 lbs. (6.5 kg) 3.25 inches (8.4 cm) 10.5 inches (27 cm) 13.5 inches (38 cm)

On/Standby

Control, Assist/Control, SIMV, CPAP, NPPV Volume Control, Pressure Control, Pressure Support, Spontaneous 0-80 breaths per minute 50 – 2000 ml 0.3 - 9.9 seconds (100 lpm) 160 lpm $1 - 99 \text{ cmH}_{2}0$ Off, 1 - 60 cmH₂O Off, 1 – 9 lpm 10 lpm during exhalation $0 - 20 \text{ cmH}_{2}0$ Easy or Hard unlock methods 1 x current settings

10-60 seconds $5 - 100 \text{ cmH}_20$ Off, 1-60 cmH₂O 60 - 85 dBA at 1 meter Red LED, audible alarm

Off, 0.1 – 99 liters 60 seconds Off, 5 – 80 bpm Off, 3 - 40 cmH₂O

 $0 - 120 \text{ cmH}_{2}0$

 $0 - 99 \text{ cmH}_{2}0$

 $0 - 99 \text{ cmH}_{2}0$

-10 - 108 cmH₂O

0 – 250 bpm

0 – 4000 ml

99:1 - 1:99

Green LED

10 - 100 lpm

0-99.9 liters

With Boot:

Power Indicators: External Power:

Green, Amber LEDs Green, Amber, Red LEDs Green, Amber, Red LEDs

Optional Accessories:

Charge Status:

Battery Level:

LTM Graphics Monitor Package - P/N 17650-001 Universal Power Supply - P/N 14546-001 Floor Stand - P/N 10611 Utility Basket - P/N 11450 Long Crossbar – P/N 11452 Circuit Support Arm – P/N 11501 Humidifier Pole – P/N 11453 Oxygen Canister Ring Set - P/N 11455 UPS Mounting Bracket - P/N 15009-001

Extended Features:

Apnea Back-up Ventilation Pressure Control and Pressure Support Variable Rise Time Pressure Control and Pressure Support Variable Flow **Termination Criteria** Pressure Support Variable Time Termination Leak Compensation

Power Specifications:

11 – 15 VDC 90-250 VAC 47 – 63 Hz

Pneumatic Specifications:

Low Pressure 02 Source:

Environmental Specifications:

Operating Temperature: Storage Temperature: **Operating Humidity:** Storage Humidity: Shock Compliance: Vibration Compliance:

5°-40°C (40°-104°F) -20°-60°C (-4°-140°F) 15% - 95% relative 10% - 95% relative MIL-STD-810E MIL-STD-810E

<80 lpm, <10 PSIG

Standards and Regulatory Compliance:

ETC, cETL Listed ISO 13485 Registration IEC 601-2-12 BSI Registered FM 75092

Pulmonetic Systems



17400 Medina Rd., Suite 100 • Minneapolis, MN 55447, USA Telephone: (763) 398-8300 • Fax: (763) 398-8400 www.pulmonetic.com

